

DECG Reusable Legplate Adapter Cable

The Connection Between Mother and Child

The DECG Reusable Legplate Adapter Cable, hereafter called the “Legplate” provides a signal path between the Fetal Spiral Electrode (FSE), Attachment Electrode, and Fetal Monitor when monitoring fetal heart rate during labor and delivery procedures.

Use a legplate that is compatible with your Fetal Monitor. Connect only compatible FSE to the Legplate.

The Legplate is compatible with:

- Philips Avalon Monitor

Article Number

AD-UFT-001	DECG Reusable Legplate Adapter Cable
------------	--------------------------------------



Follow the link for more information:
<https://www.neoventa.com/>

