

Mounting Solutions

Compatible with Stan S41 Maternal and Fetal Monitor

There are two mounting options for the S41 monitor: on a trolley or on a wall arm. Both options offer a flexible, practical, and ergonomic mounting solution. All wall arms and trolleys include an adaptor plate compatible with the S41 monitor and enables the quick-release function.

Pick and mix among trolleys, wall arms, and accessories to build a mounting system tailored entirely after your needs.

Trolley

The trolley enables easy access to necessary equipment in different locations when needed, and meets various ergonomic needs with adjustable height and position options.

- Available in fixed or variable height
- Quick release function
- Enables free positioning of the monitor
- Stable with quiet-running, easy-rolling lockable wheels and anti-static wire
- Friction-resistant tilt functionality
- A variety of accessories available



Trolley Fixed Height

Trolley with fixed height and quick-release function for the S41 monitor.

- Height: 128 cm
- Maximum storage weight: 4,5 kg
- Color: Cream white

NVT-0002-60

- 1 Trolley fixed height
- 1 Bin
- 1 Handle
- 2 Cable management clips
- 1 Anti-static wire



Trolley Variable Height

Trolley with variable height (gas-spring assisted) and quick-release function for the S41 monitor.

- Height range: 120-156 cm
- Maximum storage weight: 6,8 kg
- Color: Cream white

SBS120430

- 1 Trolley variable height
- 1 Bin
- 1 Handle
- 2 Cable management clips
- 1 Anti-static wire



Wall Arm

The wall arm offers increased free floor space, and stores accessories in a space-preserving and practical way for different ergonomic needs and room layouts.

- Available in Standard and Premium
- Quick release function
- Lateral and vertical movement functionality and tilt function
- Available in long and short wall channel lengths
- A variety of accessories available

Standard



Premium

Standard Wall Arm – Short Channel

Fixed height wall arm with lateral movement and tilt function. Includes quick-release function for S41.

- Maximum extension from wall to mounting plate: 39 cm
- Wall channel length: 48,3 cm

SBS120415/1

1 Wall arm lateral movement
1 Wall channel 48,3 cm




Standard Wall Arm – Long Channel

Fixed height wall arm with lateral movement and tilt function. Includes quick-release function for S41.

- Maximum extension from wall to mounting plate: 39 cm
- Wall channel length: 63,5 cm

SBS120415/2

1 Wall arm lateral movement
1 Wall channel 63,5 cm




Premium Wall Arm – Short Channel

Wall arm with lateral and vertical movement (gas-spring assisted) and tilt function. Includes quick-release function for S41.

- Maximum extension from wall to mounting plate: 66 cm
- Wall channel length: 48,3 cm

SBS120420/1

1 Wall arm lateral and vertical movement
1 Wall channel 48,3 cm



Premium Wall Arm – Long Channel

Wall arm with lateral and vertical movement (gas-spring assisted) and tilt function. Includes quick-release function for S41.

- Maximum extension from wall to mounting plate: 66 cm
- Wall channel length: 63,5 cm

SBS120420/2

1 Wall arm lateral and vertical movement
1 Wall channel 63,5 cm



Accessories

Trolley: ● Wall arm: ●

Basket

Storage for monitoring accessories such as cables, blood pressure cuffs.

Dimensions: W:32 cm D:19 cm H:15 cm

● RS-0001-24

● WM-0001-18B



Cord wrap

Cord management to be mounted horizontally or vertically.

● RS-0013-03E



Chart Bin

Chart bin with quick release adapter.

Dimension: W:36,2 cm. D:8,3 cm, H:22,9 cm

● WM-0001-23

● WM-0001-21



Bin

 **WM-0001-22**

Bin with removable dividers. Quick release adapter for easy removal. Included in the trolley standard package.

Dimensions: H:7,6-12,7 cm D:30,5 cm W:30 cm

 **WM-0001-27**



Quad 4-hook

 **UT-0001-11**

Hook for cable management.

Maximum load per hook: 1,27 kg each.



Follow link to find out more.

www.neoventa.com

